cefar[®] peristim^{pro}

	USERS GUIDE	3
	MANUAL DEL USUARIO	25
	GEBRAUCHSANWEISUNG	47
	MODE D'EMPLOI	69
	ISTRUZIONI PER L'USO	91
	GEbruikersHANDLEIDING	113
	BRUKSANVISNING	135

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Urge incontinence

Urge incontinence means a sudden, strong urge to urinate followed by an immediate bladder contraction, resulting in an involuntary leakage of urine. Both men and women can be affected, particularly the elderly. One reason for this condition can be a disruption in the part of the nervous system that controls the bladder.

Mixed incontinence

Mixed incontinence is a combination of stress incontinence and urge incontinence.

Faecal incontinence

Faecal incontinence, also called anal or bowel incontinence, is the impaired ability to control passage of gas or stool. There are many possible causes of faecal incontinence, the most common is injury to the anal sphincter (ring-like muscle), for instance during childbirth or surgery, or damage to the nerves that control the anal sphincters. The condition usually becomes worse as people age.

INCONTINENCE TREATMENT

Electrical stimulation via the pelvic nerves is a recognized treatment alternative for urinary incontinence. It is also proposed as a treatment method for faecal incontinence due to pelvic floor dysfunction or a poorly functioning anal sphincter.

When treating stress incontinence, the aim of the electrical stimulation is to mirror a voluntary muscle contraction and improve the function of the pelvic floor muscles. For urge incontinence, the aim is to inhibit involuntary bladder contractions by stimulating the nerves in the pelvic floor. When treating mixed incontinence, a stimulation appropriate for both urge incontinence and stress incontinence is used. For faecal incontinence, the aim is to improve bowel control by strengthening and toning pelvic floor muscles.

3. PRECAUTIONARY MEASURES

- Use only CEFAR-approved vaginal and anal probes intended for electrical stimulation.
- The stimulator may only be operated after instructions from medically trained personnel.
- The treatment is a personal ordination. Do not lend probes or stimulator to anyone else
- Only Cefar accessories should be used with the stimulator.

WARNINGS

For internal use with a vaginal/anal probe

- People with implanted electronic equipment, such as pacemakers and intracardiac defibrillators, must not be treated with CEFAR PERISTIM PRO.
- Pregnant women must not be treated with the urology programs.
- People with extra-urethral incontinence (fistula, ectopic ureter) must not be treated with the urology programs.
- People with overflow incontinence due to outflow obstacle must not be treated with the urology programs.
- People with serious retention of urine in the upper urinary tract must not be treated with the urology programs.
- People with complete peripheral denervation of the pelvic floor must not be treated with the urology programs.
- Stimulation should not take place while the user is connected to high-frequency surgical equipment. It may cause burn injuries on the tissue, as well as problems with the stimulator.
- Do not use the stimulator in the vicinity of shortwave or microwave therapy equipment, since this may affect the output power of the stimulator.
- Keep the stimulator out of reach of children.

CAUTION

For internal use with a vaginal/anal probe

- Patients with total/subtotal prolapsed uterus/vagina should be stimulated with greatest caution.
- Patients with urinary tract infections must be treated and clear of infection before starting therapy with the urology programs. Consult your doctor.
- If tissue irritation should occur, treatment should be temporarily discontinued. If problems continue, contact your health care provider. Hypersensitivity can occur in isolated cases. The problem usually disappears when probes or gel are changed to another type.
- Observe caution when using electrotherapy at the same time as the patient is connected to monitoring equipment with body worn electrodes. The stimulation might interfere with the signals to the monitoring equipment.
- Never open the battery cover during stimulation in order to avoid electrical shock.
- Turn off the stimulation or make sure that the amplitude for each channel is 0,0 mA before touching or removing the probes. Getting electrical stimulation through the fingers is unpleasant but not harmful.
- Observe caution when stimulating in the immediate vicinity of cellular phones that are switched on, since this may affect the output power of the stimulator.

For external use with surface electrodes

3.1. CONTRAINDICATIONS

- Implanted electronic devices. Do not use the device on patients who have a cardiac pacemaker, implanted defibrillator, or other implanted electronic device, because this may cause electric shock, burns, electrical interference, or death.
- TENS for Undiagnosed Pain. Do not use the device as a TENS device on patients whose pain syndromes are undiagnosed.

3.2. WARNINGS

- Consult with physician. Consult with the patient's physician before using the device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- Skin condition. Apply stimulation only to normal, intact, clean, healthy skin.
- Long term effects. The long-term effects of chronic electrical stimulation are unknown.
- Stimulation location. Stimulation over Neck or Mouth. Do not apply stimulation over the patient's neck (especially the carotid sinus) or the patient's mouth, because this could cause severe muscle spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- Stimulation across Chest. Do not apply stimulation across the patient's chest, because the introduction of electrical current into the chest may cause rhythm disturbances to the patient's heart, which could be lethal.
- Across the head. Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head, and electrodes should not be placed on opposite sides of the head.
- Stimulation over compromised skin. Do not apply stimulation over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins).
- Stimulation near cancerous lesions. Do not apply stimulation over, or in proximity to, cancerous lesions.
- Stimulation over eyes. Do not apply stimulation directly on the eyes.
- Environment. Electronic monitoring equipment. Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Bath or Shower. Do not apply stimulation when the patient is in the bath or shower. Do not apply stimulation in humid atmosphere exceeding 75% of relative humidity. Sleeping. Do not apply stimulation while the patient is sleeping.
- Driving or operating machinery. Do not apply stimulation while the patient is driving, operating machinery, or during any activity in which electrical stimulation or involuntary muscle contraction can put the patient at risk of injury.
- Electrosurgical equipment or defibrillators. Disconnect the stimulation electrodes before using electrosurgical equipment or defibrillators. Otherwise skin burns may be caused below the electrodes and the device might be destroyed.
- Magnetic Resonance Imaging. Do not wear electrode or the device during Magnetic Resonance Imaging (MRI) scans as this may result in metal overheating and causing skin burns in the area of the electrode.
- Flammable or explosive environment. Do not use the device in areas where there is a risk of fire or explosion, such as oxygen-rich environments, in the vicinity of flammable anaesthetics, etc.
- Power supply. Never connect stimulation cables to an external power supply as there is a risk of electric shock.
- Near other equipment. Do not use the device beside or stacked on top of any other equipment. If you must use it side by side or on top of another system, you should check that the device works properly in the chosen configuration.
- Miscellaneous. Electrodes for Single Patient. Do not share electrodes with other persons. All users should have individual set of electrodes to prevent undesirable skin reactions or disease transmission.
- Accessories. Use this device only with the leads, electrodes, and accessories recommended by the manufacturer. Use of other accessories may adversely affect the performance of the device or may result in stronger electromagnetic emissions or reduce the electromagnetic immunity of the device.
- No Modification. No modification of the equipment is allowed.

3.3. PRECAUTIONS

- Supervision. Use this device only under the continued supervision of a licensed practitioner. Electrode placement and stimulation settings should be based on the guidance of the prescribing practitioner.

- **Manufacturer.** The manufacturer does not take any responsibility for any electrode placements other than recommended.
- **Pregnancy.** The safety of electrical stimulation during pregnancy has not been established.
- **Skin irritation.** Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (gel). The irritation may be reduced by using an alternate conductive medium or alternate electrode placement. Some patients may experience redness under the electrodes after a session. This redness usually disappears within a few hours. Advise the patient to consult the clinician if the skin redness does not disappear after a few hours. Do not start another stimulation session in the same area if the redness is still visible. Don't scratch the redness area.
- **Heart Disease.** Patients with suspected or diagnosed heart disease should follow precautions recommended by their physicians.
- **Epilepsy.** Patients with suspected or diagnosed epilepsy should follow precautions recommended by their physicians.
- **Internal Bleeding.** Use caution when the patient has a tendency to bleed internally, such as following an injury or fracture.
- **After Surgery.** Use caution following recent surgical procedures when stimulation may disrupt the patient's healing process.
- **Over uterus.** Do not place electrodes directly over the uterus or connect pairs of electrodes across the abdomen if you are pregnant. The reason is that, theoretically, the current could affect the foetus's heart (although there are no reports of it being harmful).
- **Lack of sensation.** Use caution if stimulation is applied over areas of skin that lack normal sensation. Don't apply stimulation on patient unable to express themselves.
- **Stimulation.** The stimulator should only be used with skin electrodes intended for nerve and muscle stimulation. Muscle soreness may occur after stimulation but usually disappears within a week.
- **Hot casing or batteries.** Under extreme use conditions, some parts of the casing might reach up to 109 °F (43 °C). Use caution when manipulating the batteries right after device use or when holding the device. There is no particular health risk associated with this temperature besides your comfort.
- **Children.** Keep this device out of the reach of children.
- **Electrode Size.** Do not use electrodes with an active area less than 16 cm², as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2mA/ cm².
- **Strangulation.** Do not wrap leadwires around your neck, and keep them out of the reach of children. Strangulation may result from entanglement in the leadwires.
- **Tripping.** Care should be used to avoid tripping on lead wires.
- **Damaged Device or Accessories.** Never use the device or any of its accessories if it is damaged (case, cables, etc.) or if the battery compartment is open as there is a risk of electric shock. Carefully inspect the lead wires and connectors prior to each use.
- **Inspect Electrodes.** Inspect electrodes before each use. Replace electrodes when they begin to deteriorate or lose adhesion. Poor contact between the electrodes and the patient's skin increases the risk of skin irritation or burns. Electrodes will last longer if used and stored according to instructions on electrode packaging. Attach the electrodes in such a way that their entire surface is in contact with the skin.
- **Foreign bodies.** Do not allow any foreign bodies (soil, water, metal, etc.) to penetrate the device and the battery compartment.
- **Batteries.** Do not carry batteries in a pocket, purse, or any other place where the terminals could become short-circuited (e.g. by way of paper clip). Intense heat could be generated and injury may result. Never open the battery cover during stimulation in order to avoid electrical shock. Remove the batteries from the device if you do not

intend to use it for a prolonged period of time (more than 3 months). Extended storage of the batteries in the device might lead to batteries and device damage.

- Cable. The cable is best preserved if left attached to the stimulator between sessions. Do not jerk the cable or connection.
- Heat and Cold products. The use of heat or cold producing devices (e.g. electric heating blankets, heating pads or ice packs) may impair performance of the electrode or alter the patient's circulation/sensitivity and increase the risk of injury to the patient.
- Pulled muscles. Do not apply electrodes over pulled muscles. Using the stimulator on a previously extended muscle might further pull such muscle. The higher the stimulation intensity, the higher the risk to further overextends such muscle.
- Additional Precautions for TENS
 - TENS is not effective for pain of central origin, including headache.
 - TENS is not a substitute for pain medications and other pain management therapies.
 - TENS devices have no curative value.
 - TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- Effectiveness of TENS is highly dependent upon patient selection by a practitioner qualified in the management of pain patients.

3.4. DANGERS

- Electrodes. Any Electrode with a minimum active area of 16 cm² may be used with this device. Use of an electrode with an area less than 16 cm² can cause burns when the unit is used at higher intensities. Consult your clinician prior to using any electrode less than 16 cm².

3.5. ADVERSE REACTIONS

- Patients may experience skin irritation and burns beneath the stimulation electrodes applied to the skin.
- Patients may experience headache and other painful sensations during or following the application of electrical stimulation near the eyes and to the head and face.
- Patients should stop using the device and should consult with their physicians if they experience adverse reactions from the device.
- Precaution: Do not disconnect any stimulation cables during a session while the stimulator is switched on. Switch the stimulator off first. Always turn off the stimulator before moving or removing any electrodes during a session.
- Precaution: Do not use electrodes with an active area less than 16 cm², as there will be a risk of suffering a burn injury. Caution should always be exercised with current densities more than 2mA/cm².
- Precaution: Do not apply stimulation in the vicinity of metal. Remove jewelry, body piercings, buckles or any other removable metallic product or device in the area of stimulation. Never use the electrodes contra-laterally, i.e. do not use two pins connected to the same channel on opposite segments of the body.
- Precaution: Never carry out an initial stimulation session on a person who is standing. The first five minutes of stimulation must always be performed on a person who is sitting or lying down. In rare instances, people of a nervous disposition may experience a vasovagal reaction. This reaction is connected with fear of the muscle stimulation as well as surprise at seeing one of their muscles contract without having intentionally contracted it themselves. A vasovagal reaction causes heart to slow down and blood pressure to drop, which can make you feel weak and faint. If this does occur, stop the stimulation and lie down with the legs raised until the feeling of weakness disappears (5 to 10 minutes).
- Precaution: Sudden temperature changes can cause condensation to build up inside the stimulator. To prevent this, allow it to reach room temperature before use

- ## 2. INCREASE (left and right channel)

- Note!** Always increase the amplitude cautiously.

- ### 3. DECREASE (left and right channel)

- #### 4. PROGRAM

- ## 5. TIMER

- Initiates the timer setting.


6. PROGRAMMING/CONFIRMATION*

- Turns the stimulator into programming mode for custom programs P10–P12 when pressed for 2 seconds.
- Confirms settings in programming mode.

7. CABLE HOLDER

■■■■■

4. OVERVIEW/DISPLAY SYMBOLS

1. BROKEN CIRCUIT

Broken circuit. The reason for a broken circuit may be too high resistance or cable breakage. See chapter 10.

2. PAUSE

Paused program.

3. KEY LOCK

Activated key lock. The key lock is automatically activated if no key is pressed for 10 seconds. Deactivate the key lock by pressing the left or right **DECREASE** button.

4. PROGRAM LOCK

Activated program lock.

5. BATTERY STATUS

Empty batteries. This symbol is shown when the batteries are almost empty.

6. PROGRAM NUMBER

Selected program number.

7. REMAINING TIME

Remaining program time in minutes and seconds. Time is flashing during timer setting.

8. AMPLITUDE BARGRAPH (left and right channel)

Selected amplitude as a bargraph.

9. AMPLITUDE LEVEL (left and right channel)

Selected amplitude in mA.

10. STIMULATION/REST


Stimulation/rest indication for programs with intermittent stimulation. The upper part of the symbol is flashing during stimulation and the bottom part during rest time.

5. OPERATION

Programs 1–7 in CEFAR PERISTIM PRO should be used with the vaginal and/or the anal probe and Programs 8–9 with surface electrodes. Program 10–12 can be configured in various ways to suit individual needs and should only be used as recommended by your health care provider.


1. INSERT THE BATTERIES

Insert the batteries (see chapter REPLACEMENT OF BATTERIES).


2. PLACE THE ELECTRODES/PROBES

A1. Connect the electrodes to the cable(s).

B 1. Connect the probe(s) to the cable(s).


A 2. P8–P9: Attach the electrodes to your body


B 2. P1-P7: Insert the vaginal probe and/or the anal probe.

A 3. Connect the cable(s) to the CEFAR PERISTIM PRO.


B 3. Connect the cable(s) to the CEFAR PERISTIM PRO.

3. SWITCH THE STIMULATOR ON

Press the **ON/OFF** button, . This button can be used for terminating the stimulation at all times, even when the key lock is activated.


Always switch the stimulation off before removing electrodes from the skin.

4. SELECT A PROGRAM (P1-P12)

Select a program recommended by your health care provider by pressing the **PROGRAM** button  repeatedly until the program number is shown on the display.

Note! When selecting a program the amplitude must be set to 00.0 mA for both channels


5. START THE STIMULATION

Press the **INCREASE** button  for each channel until you reach a comfortable level of stimulation. Press and hold the button to increase the amplitude continuously.

Note! Always increase the amplitude cautiously!

Programs P1–P6: These programs have intermittent stimulation, i.e. stimulation alternates with rest time. The upper part of the stimulation/rest symbol is flashing during stimulation and the bottom part during rest time. Set the amplitude during stimulation.

The treatment time is preset, but can be changed by using the timer function, see section TIMER.

6. STOP THE STIMULATION

You can stop the stimulation before the program time has elapsed by pressing the **DECREASE** button  until the amplitudes are 00.0 mA or press the **ON/OFF** button .

An automatically activated key lock prevents unintentional changes during treatment. The key lock is activated if no button is pressed for 10 seconds. Press any of the **DECREASE** button  to deactivate the key lock.

When the program is finished, Remaining time "00:00" is flashing on the display.

7. AFTER USE

Turn the stimulation off. Remove the probe(s)/electrodes. Clean the probe(s) according to the instructions included with your probes.

The last used program is saved when the stimulator is turned off and then pre-selected the next time the stimulator is turned on.

The stimulator turns off after five minutes of inactivity to spare batteries.

To set the timer:

1. Select a program and start the stimulation, see chapter OPERATION.
2. Press the **TIMER** button  to initiate the timer function.
3. Increase the time by pressing the right **INCREASE** button . The time is increased by 1 minute each time you press the button. Decrease the time by pressing the right **DECREASE** button . The time is decreased by 1 minute each time you press the button. The time is flashing on the display, while setting the timer.
4. Press the **TIMER** button  to confirm the timer setting.

PROGRAM PAUSE

You can pause the programs up to five minutes.

To pause a program:

1. If key lock is on, press any of the **DECREASE** buttons  to deactivate it.
2. Press the **PROGRAM** button  to pause the program.

Start stimulation again by pressing the **PROGRAM** button .

If stimulation is paused for more than five minutes, the stimulator turns off automatically to spare batteries.


6. REPLACEMENT OF BATTERIES

A battery symbol is shown on the display when the batteries are almost empty . As long as the stimulator is working normally you can continue the treatment. When stimulation feels weaker than usual or the stimulator turns off, it is time to replace the batteries.

If the stimulator is not used for some time (approximately 3 months), the batteries should be removed from the stimulator.

The stimulator operates on either two non-rechargeable 1.5 V AA batteries or two rechargeable 1.2 V AA batteries, recharged in a separate charger.

REPLACING THE BATTERIES

1. Turn off the stimulator.
2. Locate the battery compartment on the back of the stimulator.
3. Remove the battery cover by sliding it downwards about 1 cm and then lifting it upwards.
4. Remove the batteries.
5. Insert the new batteries correctly according to the polarity markings (+ and -) inside the battery compartment, see picture.
6. Replace the battery cover.
7. Dispose of the exhausted batteries in accordance with local and national regulations.

NOTE! Non-rechargeable batteries may not be charged due to the risk of explosion.

7. PROGRAMS

PRESET PROGRAMS

CEFAR PERISTIM PRO offers eight preset programs for stress, urge, mixed and faecal incontinence treatment and one preset program for relief of bladder-related pain.

Programs 1-7 should be used with a vaginal and/or anal probe and programs 8-9 with surface electrodes.

In programs 1–6 intermittent stimulation is used, which means that stimulation time alternates with rest time without stimulation: stimulation/rest/stimulation. The rest time allows the muscles to recover.

Programs 1–3

STRESS INCONTINENCE

Programnumber	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P1	50 Hz	200 μ s	3 s	6 s	30 minutes, 3–5 times a week
P2	50 Hz	200 μ s	5 s	10 s	30 minutes, 3–5 times a week
P3	50 Hz	200 μ s	10 s	20 s	30 minutes, 3–5 times a week

One cause of stress incontinence can be poorly functioning perineal muscles. Electrical stimulation of the perineal muscles with a vaginal and/or anal probe makes the muscles contract. The contractions help to identify and strengthen the weakened muscles. The stimulation must be as strong as possible without being painful, and it is helpful if the patient participates actively in the muscle contractions. The intensity should be sufficient enough to create a reflexive contraction of the anus. It is suggested that this training will be combined with the patient's own exercises for pelvic floor muscles.

Programs 4–6

MIXED INCONTINENCE

Programnumber	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P4	20 Hz	200 μ s	3 s	6 s	30 minutes, 3–5 times a week
P5	20 Hz	200 μ s	5 s	10 s	30 minutes, 3–5 times a week
P6	20 Hz	200 μ s	10 s	20 s	30 minutes, 3–5 times a week

Program 4-6 is used with a vaginal and/or anal probe for treating mixed incontinence. The stimulation should be as strong as possible without being painful. It is helpful if the patient participates actively in the muscle contractions.

URGE INCONTINENCE

Programnummer	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P7	10 Hz	180 μ s	Continuous	–	30 minutes, 2–5 times a week

Normally the bladder muscle is relaxed between evacuations. In certain kinds of urge incontinence, however, the bladder is unstable and sudden uncontrolled bladder contractions occur during the "resting phase". Urge incontinence can be treated with low-frequency constant stimulation, which has a relaxing effect on the overactive bladder. Program 7 is intended for urge incontinence treatment and should be used with a vaginal and/or anal probe. The stimulation should be as strong as possible without being painful.

FAECAL INCONTINENCE

Programnumber	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P1	50 Hz	200 µs	3 s	6 s	30 minutes, 3–5 times a week
P2	50 Hz	200 µs	5 s	10 s	30 minutes, 3–5 times a week
P3	50 Hz	200 µs	10 s	20 s	30 minutes, 3–5 times a week
P4	20 Hz	200 µs	3 s	6 s	30 minutes, 3–5 times a week
P5	20 Hz	200 µs	5 s	10 s	30 minutes, 3–5 times a week
P6	20 Hz	200 µs	10 s	20 s	30 minutes, 3–5 times a week


Faecal incontinence can be a result of e.g. weakened pelvic floor muscles, weakened or poorly functioning anal sphincter muscles or damage to the nerves controlling the anal sphincter muscles. In these cases electrical stimulation can be used with an anal probe to increase the tonus of the pelvic floor and the anal sphincter. The stimulation should be as strong as possible without being painful. The intensity should be sufficient enough to create a reflexive contraction of the anus. It is helpful if the patient participates actively in the muscle contractions.

TENS (HIGH FREQUENCY) – PAIN RELIEF

Programnummer	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P8	80 Hz	180 μ s	Continuous	–	At least 30 minutes daily

High-frequency TENS stimulation is an effective treatment for relief of both acute and chronic pain. External stimulation with surface electrodes can be used for diagnosed pain conditions in the lower part of the abdomen, as in cystitis.

The electrodes should be placed in the painful area over the pubic region, see electrode placement picture.

Program 9


TENS (BURST) – URINARY INCONTINENCE

Programnumber	Frequency	Pulse width	Stimulation time	Rest time	Recommended treatment time
P9 Sacral nerves stimulation	2 Hz	180 µs	Continuous	–	At least 30 minutes daily
P9 Acupuncture pointsstimulation	2 Hz	180 µs	Continuous	–	30 minutes per session during at least one month. Start with 2–3 times a week and then reduce the treatment frequency.

Electrical stimulation for urinary incontinence is usually applied by a probe inserted in the vagina and/or anus. An alternative method is to use TENS to stimulate areas of the skin that are innervated by the same spinal cord segment (S2–S3) as the bladder and urethra. The stimulation is applied through surface electrodes placed on the skin between the anus and the genitals, alternatively in the lumbar region. In treatment of incontinence where the sacral nerves are stimulated, the electrodes should be placed in the region of S2–S3, see electrode placement picture. The stimulation should be strong enough to create a reflexive contraction of the anus.

Electrical stimulation of the posterior tibial nerve, via acupuncture points, has also proven to have a positive effect.

In treatment of incontinence where stimulation is applied over acupuncture points, the electrodes should be placed over SP6 and behind/beneath the medial malleolus, see electrode placement picture

CUSTOM PROGRAMS (P10-P12)


With the CEFAR PERISTIM PRO it is possible to create and store three custom programs for individual needs. To create a custom program, follow the programming procedure below. To use a custom program, follow the instructions in the chapter OPERATION. The program time is preset to 30 minutes but can be changed with the timer function, see section TIMER.

Programming

1. Press the **PROGRAM** button  repeatedly until Program 10, 11 or 12 is shown on the display. Select one of these programs.
2. Press the **PROGRAMMING/CONFIRMATION** button  for 2 seconds to enter the programming mode.
3. Press the **INCREASE (OR DECREASE)** button   to see the two different stimulation types alternate on the display: C, I. The programming chart below shows the setting possibilities for each stimulation type.

- A custom program can be changed by redoing the programming.

PROGRAMMING CHART

To activate/deactivate the program lock:

- ## COMPLIANCE

1. Turn the stimulator ON.
2. Press the **TIMER** button  and the right **DECREASE** button  simultaneously for 2 seconds.
3. The left side of the display shows the usage time in hours and the right side in minutes. To reset the usage time, press the right **DECREASE** button  for 2 seconds.
4. Wait for 5 seconds or press **TIMER** button .
5. The left side of the display shows the total usage time in hours and the right side in weeks. The total usage time cannot be reset.
6. Wait for 5 seconds or press the **TIMER** button  to exit compliance mode.



8. ACCESSORIES

You can use the CEFAR PERISTIM PRO together with vaginal and/or anal probe. The probes can be either autoclavable or non-autoclavable. Always read and follow the care and maintenance instructions delivered with the probes.

If you use the CEFAR PERISTIM PRO with surface electrodes the electrodes will eventually wear out and need to be replaced. It is recommended to replace the electrodes after approximately 20–40 times of usage.

CEFAR PERISTIM PRO comes with a neck strap and a belt clip, allowing you to wear the stimulator around your neck or on your belt and have your hands free during treatment.

For purchase information, contact your Cefar dealer or visit www.cefarcampex.com.

9. CARE INSTRUCTIONS

Taking care of and cleaning the Cefar-Compex equipment is simple with the following instructions:

- For care and cleaning instructions of probes, see probe package.
- Self-adhesive multi-use electrodes are re-moistened if necessary with a few drops of water and kept air-tight (in a plastic bag) on protective paper when they are not in use.
- Rinse the carbon rubber electrodes and the skin with water after use. Do not use detergent for the electrodes.
- Never expose the stimulator to water. Wipe it off with a damp cloth if necessary.
- Do not jerk cables or connections.
- The cables are best preserved if left attached to the stimulator between sessions.
- No other maintenance of the device is required. Service life of the device can vary depending on usage conditions. Typical service life is 7 years.
- Never service the device while in use
- The device should be operated in temperatures between 10°C and 40°C, atmospheric pressures between 50 and 106 kPa, and relative humidity between 30% and 75%.
- The device should be transported and stored in temperatures between -40 °C and 70 °C, atmospheric pressures between 50 and 106 kPa and relative humidity between 10% and 90%.
- Contact manufacturer for assistance in setting up, using or maintaining the equipment or report events.

10. TROUBLESHOOTING

THE STIMULATION DOES NOT FEEL THE SAME AS USUAL

- Check that all settings are correct (chapter OPERATION) and make sure that the electrodes are correctly placed.
- **Electrodes:** Slightly change the position of the electrodes.

THE STIMULATION FEELS UNPLEASANT

- Check that the probe/electrode cable is properly connected to the stimulator.
- **Probes:** Make sure that the probe is correctly placed with good contact. In case of irritation/inflammation in the vagina or anus, please contact your health care provider.
- **Electrodes:** The skin is irritated. For advice on skin care, see chapter PRECAUTIONARY MEASURE.
- **Electrodes:** The electrodes begin to lose their stickiness and do not stick properly to the skin. Moisten the adhesive surface with a few drops of water before placing on the skin.
- **Electrodes:** The electrodes are worn out and need to be replaced.
- **Electrodes:** Slightly change the position of the electrodes.

THE STIMULATION FEELS WEAK OR NOT AT ALL

- Check if the batteries need to be replaced, see chapter REPLACEMENT OF BATTERIES.
- **Electrodes:** Electrodes are too old and need replacement.

THE BROKEN CIRCUIT SYMBOL IS SHOWN ON THE DISPLAY SYMBOL

The broken circuit symbol indicates that the resistance is too high, or that a cable is broken.

- A too high resistance can be caused by a bad connection between the electrodes/probes and your skin/tissue, or that the electrodes need to be replaced.
- A cable breakage can be checked by pressing the cable's pins against one another while increasing the amplitude for the corresponding channel to 11 mA. If the amplitude now drops to 00.0 mA and  starts flashing, the cable needs to be replaced.

Note! Never increase the amplitude above 20 mA when you check for cable breaks, since this can damage the stimulator.

THE STIMULATOR IS NOT WORKING

If the error symbol appears on the display when you start the stimulator, it means that the stimulator is broken and needs to be repaired.

Note! Do not use the stimulator – contact your Cefar dealer.

Cefar-Complex will only be responsible for service and repairs performed by Cefar or a distributor appointed by Cefar-Complex.

KEY TO THE SYMBOLS

Reference number, part number

Follow instruction for uses

Type BF applied parts

Council Directive 2012/19/EU concerning Waste Electrical and Electronic Equipment (WEEE) requires not to dispose of WEEE as municipal waste. Contact your local distributor for information of the unit and accessories

Keep the device dry

Minimum and maximum temperature indications to respect

This equipment complies with all requirements of the Medical Device Directive (93/42/EEC)

Manufacturer name and address, manufacturing year

Power/Pause

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CEFAR PERISTIM PRO Electrotherapy System is intended for use in the electromagnetic environment specified below. The customer or the user of the CEFAR PERISTIM PRO Electrotherapy System should assure that it is used in such an environment

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Risk assessment on the CEFAR PERISTIM PRO Electrotherapy System indicates the compliance levels claimed are acceptable when ESD-precautionary measures are taken.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	Not Applicable - Battery powered Not Applicable - signal lines less than 3 meters	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+ 1kV differential mode (line to line) + 2kV common mode (line to ground)	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0,5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Not Applicable - Battery powered	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CEFAR PERISTIM PRO Electrotherapy System requires continued operation during power mains interruptions, it is recommended that the CEFAR PERISTIM PRO Electrotherapy System be powered from an uninterrupted power supply.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the AC mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The CEFAR PERISTIM PRO Electrotherapy System is intended for use in the electromagnetic environment defined below. The customer or the user of the CEFAR PERISTIM PRO Electrotherapy System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	[V1] V, where V1 = 3V	Portable and mobile RF communications equipment should be used no closer to any part of the CEFAR PERISTIM PRO Electrotherapy System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \frac{3.5}{f} \sqrt{P}$ 80 MHz to 800 MHz $d = \frac{3.5}{f} \sqrt{P}$ 800 MHz to 2,5 GHz $d = \frac{7}{f} \sqrt{P}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CEFAR PERISTIM PRO Electrotherapy System is used exceeds the applicable RF compliance level above, the CEFAR PERISTIM PRO Electrotherapy System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CEFAR PERISTIM PRO Electrotherapy System.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V] V/m.

Recommended separation distances between
portable and mobile RF communications equipment and the CEFAR PERISTIM PRO Electrotherapy System

The CEFAR PERISTIM PRO Electrotherapy System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CEFAR PERISTIM PRO Electrotherapy System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CEFAR PERISTIM PRO Electrotherapy System as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter P (W)	Separation distance according to frequency of transmitter d (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d = \frac{3,5}{[V]}$ √P (where V1 = 3V)	$d = \frac{3,5}{[E]}$ √P (where E1 = 3V/m)	$d = \frac{7}{[E]}$ √P (where E1 = 3V/m)
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

FCC REQUIREMENTS

Part 15 of the FCC Requirements	
This device complies with Part 15 of the FCC Rules. Operation is subject to the following 2 conditions:	• This device may not cause harmful interference • This device must accept any interference received, including the interference that may cause undesired operation.
FCC ID	T9J-RN42
Contains Transmitter Module IC	6514A-RN42

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Electrical stimulation: A physiologic approach to the treatment of urinary incontinence.
Journal of Urologic Clinics of North America; 2, Vol 18: 393-407, 1991

Walsh et al. Non-invasive antidromic neurostimulation, a simple effective method for improving bladder storage. *Neurol Urodyn* 2001;20(1):73-84

Zöllner-Nielsen M, Samuelsson S.M

Maximal electrical stimulation of patients with frequency, urgency and urge incontinence. Report of 38 cases. *Acta Obstet Gynecol Scand* 1992;71:629-631

Sand P.K, Richardson D.A Et al.

Pelvic floor electrical stimulation in the treatment of genuine stress incontinence: A multicenter, placebo-controlled trial. *Am J Obstet Gynecol* Vol 173, No 1, 1995

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Electrostimulation for urinary incontinence is often performed through an electrode that is inserted into the vagina and/or anus. An alternative method is to use TENS and stimulate the skin on the thighs with nerves from the same sacral segment (S2–S3) as the bladder and urethra (urethral). The stimulation of the thighs through skin electrodes should be placed on the skin between the anus and the genitalia, alternatively in the buttocks. In the treatment of incontinence the sacral nerves are stimulated by the electrodes placed in the area S2–S3, the electrode placement on the skin. The intensity should be a reasonably high for to create a reflex massaging effect on the anal sphincter.

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© 03/2020 DJO - 1566 - Rev C